

Paired-Tag: Single-Cell Epigenetic Profiling for Biomarker Discovery

Christopher Hartl, Haruhiko Ishii, Camila Saldanha, Pei Lin

Epigenome Technologies

Epigenome
Technologies



Epigenetics holds key biomarkers, but has been challenging to study

Epigenome as Biomarker

- Epigenetic dynamics, such as DNA methylation, histone modifications, and transcription factor binding are essential biomarkers for understanding disease mechanisms and cellular states due to their stable and regulatory roles in gene expression.
- Traditional methods, such as sorting, suffer from biases and lack the resolution needed to detect, cell-specific epigenetic variations.
- These limitations hinder the accurate identification and validation of reliable biomarkers, impeding advancements in diagnostics and personalized medicine.
- Addressing these challenges is crucial for harnessing the full potential of epigenetic biomarkers in clinical and research settings.





Figure: Epigenetics plays a central role in cell biology, mediating the long-term effects that environment and aging have on cellular metabolism and function – which in turn feed back to impact the epigenome.

Paired-Tag Enables Cell Resolution, Holistic Epigenetic and Transcriptomic Profiling



The Paired-Tag technology combines single-cell transcriptomics with scCut&Tag (scRNA + scCut&Tag) providing the ability to obtain cell-specific epigenetic biomarkers from heterogeneous cell populations. Paired-Tag can map histone modifications, chromatin remodelers, and transcription factor binding, for a complete picture of cellular epigenetic state. It is compatible with several single-cell platforms, requires no special processing, and can be multiplexed on-platform for higher experimental throughput.




Figure: UMAPP/LSI embedding of 74,770 PBMCs. Paired-Tag (RNA component)

Paired-Tag allows us to capture cell-specific profiles from heterogeneous tissues without sorting, and to identify cell types, as it easily integrates with existing single-cell atlases.



Figure: Single-cell Paired-Tag of human PBMC. Paired-Tag was applied to PBMC from one donor and profiled for the activation mark H3K27ac. Approximately 12,000 cells were recovered.

Epigenetic analysis: H3K27ac Signac scCut&Tag profiles merged across donors

Cluster identification: LSI identifies 8 clusters, including an additional T-cell cluster not seen in RNA data

Discovery of rare subtype: Rare epigenetically-defined T-cell subtype (~5% overall) was not detectable by RNA alone

Significance: Demonstrates power of single-cell epigenetic profiling to discover novel biomarkers

PBMC Profiling: First Chromatin State Atlas of Immune Cells

The immune system must respond to and remember environmental exposures and signals. Responses are driven by highly dynamic epigenetic mechanisms that have not been examined at single-cell resolution.

We profiled PBMCs from 4 donors for 4 chromatin marks using Cut&Tag and Paired-Tag, recovering >70,000 single-cell profiles.




Figure: UMAPP/LSI embedding of 74,770 PBMCs. Paired-Tag (RNA component)




Figure: Pseudobulk and single-cell H3K27me3 and H3K27me3 tracks at major cell type markers

Droplet Paired-Tag integrates with existing single-cell atlases




Figure: Integrative embedding of 4 Paired-Tag experiments (RNA only) with the SEA-AD single cell consortium data. Though multiple histone marks are targeted, only the RNA component is needed for seamless data integration.


Paired-tag makes it easy to annotate cell types and to integrate cell-level epigenetics with existing single-cell studies.

Cell annotation is difficult even for single-cell ATAC data, and presents a significant challenge for single-cell Cut&Tag.

By leveraging transcriptomics, Paired-Tag data can be integrated with existing cell atlases, streamlining analytics and maximizing data value.

Highly-reproducible scCut&Tag signatures even for TFs and PRC

Genome-wide profile correlations in C4-2 cells between bulk and single-cell Cut&Tag for a histone mark, a transcription factor, and a chromatin remodeler, demonstrating that single-cell profiles are strongly similar to bulk across each target type.



Consistency of Paired-Tag was ensured by comparing genome-wide profile correlations in C4-2 cells between bulk and single-cell Cut&Tag across H3K27me3 (histone), AR (Androgen Receptor) (TF), and EZH2 (chromatin).

Single-cell profiles closely mirror bulk for each target, highlighting strong reproducibility and consistency.

Paired-Tag retains the high accuracy of Cut&Tag down to the single-cell level, providing reliable epigenetic profiling of diverse chromatin targets across cell types.

High-efficiency single-cell profiles

To assess the profiling efficiency of Paired-Tag, the number of fragments (UMI) per cell and per modality were compared. Paired-Tag was performed on human cortical samples, targeting H3K27ac, H3K27me3 and EZH2.

Across all targets, 1000s of UMI/cell were recovered, sufficient to perform LSI-based (Signac) cell clustering on Paired-Tag data, and can produce the coverage equivalent of bulk Cut&Tag with as few as 5,000 cells.




Figure: Recovered DNA and RNA fragments from three Paired-Tag experiments. Colored areas: Ideal (red) >1250 RNA UMI, >2000 DNA UMI; Good (yellow) >1000RNA >1250 DNA; RNA-Low (Green), DNA-Low (Blue), and Empty (purple)

Simple workflow can be multiplexed for up to 4x sample throughput




Figure (top): The Paired-Tag Workflow
Figure (bottom): The Paired-Tag Workflow

Easy-to-use: Short workflow, existing platform compatibility
Flexible: Off-the-shelf antibodies tailor your experiment to any target
High-throughput: Multiplexing enables up to 4 simultaneous experiment per 10X lane
Powerful: Demonstrated on TFs, chromatin remodelers, and histones

Towards chromatin-based epigenetic biomarker discovery

Single-cell profiling of PBMC identifies rare epigenetic T-cell subtype




Figure: UMAPP/LSI embedding of PBMC RNA from 2 donors and 10 Paired-Tag experiments.

Profiling Method: Utilized Paired-Tag to analyze PBMCs from 2 donors, targeting 4 epigenetic marks

Data Processing: Cell Ranger, Seurat, BigSur; annotated using CellTypist

Identified Cell Types: T cells (2 naive, CD8, NK, MAIT), B cells (naive + memory), Monocytes (classical, non-classical), Dendritic (pDC)


Epigenetic analysis: H3K27ac Signac scCut&Tag profiles merged across donors

Cluster identification: LSI identifies 8 clusters, including an additional T-cell cluster not seen in RNA data

Discovery of rare subtype: Rare epigenetically-defined T-cell subtype (~5% overall) was not detectable by RNA alone

Significance: Demonstrates power of single-cell epigenetic profiling to discover novel biomarkers

Neuronal Epigenetic De-repression is a key epigenetic marker of Alzheimer's



20 (10 AD, 20 CTL) brains profiled targeting 4 major marks using SplitPool Paired-Tag, yielding ~200K nuclei (50K/mark).

Repressed regions were defined on the basis of H3K9me3 deposition in healthy controls and external ChIP-seq profiling of brain and neuronal cell lines.

In excitatory neurons H3K27ac deposition was observed to a far greater extent in AD samples than healthy controls.

Epigenetic de-repression has been posited as a mechanism of AD pathology. This work highlights the power of single-cell epigenetics in uncovering novel biomarkers.




Figure: Average target occupancy (CPM) at repressed regions (x-axis) across various cell lines. (y-axis) The difference in CPM (AD - CTL) for the corresponding element. Notably, the more overall H3K27ac deposition in excitatory neurons, the larger the AD/CTL difference, highlighting the AD-specificity of the epigenetic state.

Conclusions

- Joint Profiling of Single-Cell Cut&Tag and RNA:** Provides a detailed view of DNA-binding events, similar to ChIP-seq, but with single-cell resolution.
- Incorporates RNA Profiling:** Identifies cell types and maps interactions between enhancers and their target genes.
- Valuable Insights into Epigenetic Drivers:** Uncovers key cellular behaviors and processes, spanning contexts from development to disease.
- Future Potential:** Innovations like Paired-Tag pave the way for deeper biological discoveries and promising therapeutic breakthroughs.

Supported by:
• 1R41MH128993
• 1R43AG079691
• 2R44GM146330
• 2R44AG079691

Epigenome
Technologies

